



Molly[®]
Modular Platform
Technology



Volume 

Viscosity 

Device Type

Disposable autoinjector

Primary Container

1.0 mL long pre-filled sryinge with staked needle

Dosage Type

Fixed | Full contents

Route of Administration

Subcutaneous

Activation Method

Needle cover

Injection Feedback

Audible click at beginning and toward end of injection, viewing window

Handling Instructions



Remove cap.



Press and hold against injection site to activate.



Discard in sharps container.

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications