



Molly[®]
Modular Platform
Technology



Volume 

Viscosity 

Device Type

Disposable autoinjector

Primary Container

2.25 mL pre-filled syringe with staked needle

Dosage Type

Fixed | Full contents

Route of Administration

Subcutaneous

Activation Method

Needle cover

Injection Feedback

Audible click at beginning and toward end of injection, viewing window

Handling Instructions



1
Remove cap.



2
Press and hold against injection site to activate.



+

Discard in sharps container.

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic cap removal
- ✓ Robust anti-roll design



The SHL Medical logo and Molly are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, and other jurisdictions. Final product specifications are subject to change based on customer requirements. Contact info@shl-medical.com for more information.

www.shl-medical.com