









Device Tupe

Primary Container Dosage Type

Activation Method Injection Feedback Disposable autoinjector

3.0 mL standard ISO cartridge with pre-installed needle (NIT®)

Fixed Dose | Optional automated priming and titration

Route of Administration Subcutaneous I Intramuscular

Needle cover

Continuous audible clicks during injection, viewing window

* For suspension drugs that require a priming step - hold device with cap upright during cap removal.

Copyright © 2024 SHL Medical AG. All rights reserved. Version 2.6 Information accurate at time of print and is subject to change without notice. Contact info@shl-medical.com for more information.

Handling Instructions



Unscrew cap until needle cover automatically extends*. Then remove cap.



Press and hold against injection site to activate.



Discard in sharps container.

Key Features

- √ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Continuous audible and visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle

- ✓ Compact design
- ✓ Customizable needle length and gauge
- ✓ Passive locking needle cover
- ✓ Customized to client specifications

